



Manual Treadmill 939908/09/10

Cardio Pulmonary



Lode BV
Zernikepark 16
NL 9747 AN
Groningen
The Netherlands
t: +31 50 571 28 11
f: +31 50 571 67 46
@: ask@lode.nl
www.lode.nl

Table of contents

1 - Intended Use	3
2 - Precautions	3
3 - Contra Indication	7
4 - Validity	7
5 - Definition of parts	8
6 - Install your treadmill	8
6.1 - Emergency stop installation	9
7 - Maintenance	10
8 - Troubleshooting	10
9 - Cleaning	11
9.1 - Cleaning	11
10 - List of symbols used	11
11 - Electromagnetic Compatibility (EMC)	13
12 - Specification remarks	16
13 - Accessories	18
14 - Serviceparts	19
15 - Specifications	21



1 - Intended Use

This device is intended to be used as a stress test device in a medical environment. The main goal of the use of the devices is to create reproducible stress tests. With cycling ergometers typically workload (watt) is imposed. With treadmills typically speed (km/h) and inclination (% grade) are imposed. This product is designed both for manual operation and for control by external ECG-, pulmonary equipment. Most likely this device will be used in conjunction with another medical device to obtain other important physiological data, allowing a physician to evaluate a test subjects physical status. The product may also be used for rehabilitation or active aging therapy. As such the intended use of the device solely cannot be defined specifically.

The device has to be operated under the supervision of well-trained medical specialists in the field of use.

The ultimate judgment whether a test subject should undertake a stress test with the device must be made by the responsible medical specialist, based on the limitations of each individual, the medical history and all other applicable circumstances. Neither the manufacturer nor its distributors assume any responsibility for the final use of its equipment.

2 - Precautions



Read this manual before using the device and follow it carefully.



The operator should instruct the test subject prior to performing an exercise protocol. If, at any time during exercise, the test subject feels faint, dizzy, or experiences pain, stop the test and he or she should be consulted by the physician.



Care should be taken in mounting or dismounting the ergometer. Be aware of feet when replacing the ergometer.



The operator should not touch accessible parts and the patient simultaneously.



The test subject should always wear the belt of the lanyard safety belt.



Replacement of the power supply cord should be installed by authorised service personnel, instructions for correct connection and anchoring must be ensured, securing that the requirements of clause 8.11.3 of IEC 60601 1:2005 are met.



Set up and operate the device on a solid level surface.



The test subject should not wear loose or dangling clothing while using the device.



Service of this device is restricted to factory trained personnel only.



Do not jump onto the rotating belt. Do not jump off the rotating belt (not even to the front). Do not stop moving on the running belt. Do not turn around on the running belt. Do not move sideways or backwards. Do not make movement that could get you out of balance.



The test subject shall wear clothes when using the harness or the lanyard belt to prevent skin irritation.



Replacement of parts can only be done by the manufacturer or designated service personnel.



Read all warnings posted on the device.



Inspect the device for worn or loose components prior to use. Tighten / replace any loose or worn components prior to use.



Do not place any loose parts or objects on the device before and during use.



The equipment has a safety earth (ground) connection and must be connected to a (grounded) wall socket with protective earth to avoid the risk of electric shock. The functional earth connection is for potential equalisation only.



The operator must keep away from the footrail while lowering.



The test subject should wear proper sports or running shoes (no spikes).



If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.



Cleaning and user maintenance shall not be made by children without supervision.



Children being supervised are not to play with the appliance.



Children shall not play with the appliance.



Heart rate monitoring systems, BP measurements or SpO₂ monitoring may be inaccurate. Over exercise may result in serious injury or death. If the test person feels faint stop exercising immediately.



Keep away from the rollers.



No modification of this equipment is allowed.



Not suitable for use in the presence of flammable anaesthetics.



The appliance is not to be used by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge.



The test subject should wear proper sports clothing. The operator should check possible entrapment of clothing, fingers or feet, before the start of an exercise test.



This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.



This device should only be sold by, or under the supervision of authorized persons.



Your device may be executed with heart rate measurement, blood pressure measurement, and / or SpO₂ (blood saturation) measurement. These measurements may be inaccurate, depending on use circumstances. These measurements are indicative and cannot be used for diagnostic purposes.

3 - Contra Indication

The device is to be operated by classified personnel only. As stated in the intended use, the device is intended to be used in a medical environment. During the intended use the test subject will deliver energy. Application of the wrong dosis of energy could lead to permanent damage of the test subject health. Only use this device under supervision of a medical physician.

4 - Validity




This manual covers all OEM treadmill versions starting with the following serial number:

- 939907: S/N 20150001 -
- 939908: S/N 20150001 -
- 939909: S/N 20150001 -
- 939910: S/N 20150001 -

Date of issue: July 1, 2015

5 - Definition of parts



1. Service Connector
2. Communication Interface Module (CIM)
3. On/Off switch
4. Equipotential bonding plug
5. Safety symbols tag
6. Serial number tag
7. Motor compartment for height adjustment and drive
8. Height adjustment (optional)
9. Walking surface (APPLIED PART)
10. Front bar
11. Connector for external emergency stop
12. Magnetic contact for safety belt

6 - Install your treadmill



1. While unpacking the Lode treadmill ergometer, leave the transport straps on both ends in place.
2. Place the treadmill in the location where it will be used.
3. Remove Both transport straps and remaining packaging materials.
4. The power supply cord should be installed by authorised service personnel and instructions for correct connection and anchoring must be ensured, securing that the requirements of clause 8.11.3 of IEC 60601-1:2005 are met.
5. Remove all objects from the walking surface.
6. Switch on the unit with the on/off button and wait until the self-test is carried out. The product is ready for use after you have heard a series of “beeps”.

6.1 - Emergency stop installation



When your treadmill is delivered with an emergency stop, please note the following:

The plug of the emergency stop button needs to be placed in the back connector of the console on the front handrail.

Some emergency stop buttons have a connector for another emergency stop. You can place the spare dummy there in case you have one button.

Also make sure that the magnet of the emergency lanyard is connected to the console. Otherwise the treadmill will not operate.



7 - Maintenance

Maintenance should be carried out on a regular and planned basis. We recommend to check the unit annually. This may be done by your local dealer. It is also recommended that a record of the service history is kept for all activities relating to service and maintenance.

Maintenance and all repairs should only be carried out by an authorized agency. The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorized persons.

The check up and/or technical maintenance must be carried out conform the procedure described in the service manual of the unit.

Opening of the equipment by unauthorized agencies is not allowed and will terminate any claim to warranty.

Lode will make available on request circuit diagrams, component part lists, descriptions, calibration instructions, or other information that will assist authorized service personnel to repair those parts of the ergometer that are designated by Lode as repairable by authorized service personnel.

8 - Troubleshooting

For a complete overview of possible error message and the solutions belonging to these messages, please consult the service manual of the ergometer.



9 - Cleaning

The belt can be cleaned with a sponge, with little lukewarm water. Allow the belt to dry thoroughly after cleaning before re-using the treadmill.

9.1 - Cleaning

Before cleaning and other maintenance of any part of the device, first switch off the device and disconnect it from the mains. The surfaces can be cleaned with a damp cloth and a liquid (non-abrasive) cleaning product.

General Cleaning

Wipe the device surface down with a cloth moistened with soap water or a disinfectant. The cloth should not be dripping wet; do not allow liquids to enter the device.

Cleaning the Saddle

Clean the saddle with a soft and dry or moist cloth. Disinfectants used should not contain any alcohol.

Cleaning the Upholstery (e.g. couch ergometer)

Wipe the upholstery down with a soft cloth moistened with soap water. The cloth should only be moist and not dripping wet. If the cleaning agents and disinfectants used are caustic or contain alcohol, they may damage and/or discolor the upholstery.

If you are not certain about the discolouring effect of a cleaning agent, you may try a little on a part of the product that is not visible during normal use.

Lode cannot advise a specific cleaning agent, since local recipes may differ.

10 - List of symbols used



Read manufacturer's guide, advises and instructions and manual



FCC 18 RF ISM tested



Symbol for collection, treatment, recycling and disposal of waste electrical and electronic equipment (WEEE) as set out in Directive 2002/96/EC of 27 January 2003 of the European Parliament and of the Council on waste electrical and electronic equipment are necessary to reduce the waste management problems linked to the heavy metals concerned and the flame retardants concerned



On-off Switch



Emergency Stop



Tested and certified to U.S. and Canadian National Standards by a NRTL, viz TÜV Rheinland. Compliance for U.S. and Canadian markets.



Class I MEDICAL EQUIPMENT Type B Electrical Safety IEC 60601-1



Safety notes



Potential Equalization Conductor



**Risk by moving parts:
Rotating rollers & belt
Lowering foot rail**

Trapping zone. Warning: do not jump onto the rotating belt. Do not jump off the rotating belt (not even to the front). Do not stop moving on the running belt. Do not turn around on the running belt. Do not move sideways or backwards. Do not make movement that could get you out of balance. Keep away from rollers. The operator must keep away from

the footrail while lowering.

 RS232



 External Control Connector

 USB

11 - Electromagnetic Compatibility (EMC)



The Lode equipment is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the Lode equipment is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR11	Group 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The equipment is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations / Flicker emissions IEC 61000-3-3	Complies	
Harmonic Emissions IEC 61000-3-2	Not applicable	Equipment for professional use with a total rated power greater than 1 kW and is excluded from this requirement

The Lode equipment is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the Lode equipment is used in such an environment.


Immunity Test	IEC 60601-Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) EN 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EN 61000-4-4	Not Applicable	Not Applicable	Conducted only on ports interfacing with cables whose total length, according to the manufacturer's functional specification, may exceed 3 meters
Surge EN 61000-4-5	Not Applicable	Not Applicable	Conducted only on ports interfacing with cables whose total length, according to the manufacturer's functional specification, may exceed 3 meters
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	Not applicable	Not applicable	Equipment for professional use with a total rated power greater than 1 kW and is excluded from this requirement
Power frequency (50/60	3 A/m	3 A/m	Power frequency magnetic



Hz) magnetic field EN 61000-4-8			fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.
------------------------------------	--	--	---

NOTE: *Ut* is the AC mains voltage prior to application of the test level.

The Lode equipment is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the Lode equipment is used in such an environment.

Immunity Test	IEC 60601-Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Lode equipment including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: $d = 1,2 \sqrt{P}$</p>
Radiated RF 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	<p>$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a, should be less than the compliance level in each frequency range. b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 



NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment.
- b. Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended Separation Distances between Portable and Mobile RF Communication Equipment and the Lode equipment

The Lode equipment is intended for use in the electromagnetic environment on which radiated RF disturbances are controlled. The customer or the user of the Lode equipment can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Lode equipment as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output Power Transmitter in Watts W	Separation Distance in Meters (m) according to frequency of Transmitter		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,7	3,7	7,37
100	11,7	11,7	23,3

For transmitters rated at a maximum output power not listed above, the recommended separation distance [d] in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE:

These guidelines may not apply in all instances. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

12 - Specification remarks













The specification in this manual is of the basic version of the product at the time of the publication of this manual. The specification of all actual versions can be found on the Lode website. www.lode.nl

Changes to specifications can be made without prior notice.




Mentioned speed at the specifications in this manual is the standard speed supplied on the treadmill. When optional speed upgrade or speed downgrade is installed on the treadmill, the actual speed of your treadmill may differ from the speed mentioned at the specifications.








13 - Accessories

<p>Emergency Stop Button</p> <p>Partnumber: 945804</p> <p>Ultimate safety</p>  <p>Emergency Stop Button</p>	<p>Handrails, Side - Extension</p> <p>Partnumber: 939816</p> <p>Extended support for test subjects</p>  <p>Extension of the side handrails</p>	<p>Entrance plate</p> <p>Partnumber: 938809</p> <p>Even easier entrance to the treadmill</p>  <p>Entrance plate</p>	<p>Universal treadmill Arm Support</p> <p>Partnumber: 945805</p> <p>Comfort for both test subject and test</p>  <p>Universal treadmill Arm Support</p>	<p>Colour display 3.5" - single display</p> <p>Partnumber: 945810</p> <p>Clear feedback</p>  <p>Colour display 3.5" - single display</p>
<p>Control Unit with 7" touch screen for treadmill</p> <p>Partnumber: 945814</p> <p>Multifunctionality</p>  <p>Control Unit with Touchscreen</p>	<p>Programmable Control Unit with 7" touch screen for ...</p> <p>Partnumber: 945815</p> <p>Programmable</p>  <p>Programmable Control Unit</p>	<p>Colour Display 3.5" - 2nd screen</p> <p>Partnumber: 945819</p> <p>Multifunctionality</p>  <p>Colour Display 3.5"</p>	<p>Blood Pressure Measurement with ECG trigger for ...</p> <p>Partnumber: 945824</p> <p>with ECG trigger</p>  <p>Blood Pressure Measurement with ECG trigger</p>	<p>SpO2 for control unit with touch panel (extra long ...)</p> <p>Partnumber: 945822</p> <p>Oxygen saturation</p>  <p>SpO2 for control unit with touch panel</p>













<p>USB to Serial converter</p> <p>Partnumber: 226012</p> <p>Easy connection</p>  <p>USB to Serial converter</p>	<p>RS232 cable</p> <p>Partnumber: 930911</p> <p>Easy connection</p>  <p>RS232 cable</p>	<p>Table foot for Control Unit with Touchscreen</p> <p>Partnumber: 945818</p> <p>Control from any desired place</p>  <p>Support for Control Unit with Touchscreen</p>
---	---	---

14 - Serviceparts

<p>3,5" Control Unit front panel</p> <p>Partnumber: 945110</p> <p>Original Lode Service part</p>  <p>3,5" Control Unit front panel</p>	<p>7" Control Unit front panel</p> <p>Partnumber: F945105</p> <p>Original Lode Service part</p>  <p>7" Control Unit front panel</p>	<p>7" Programmable Control Unit front panel</p> <p>Partnumber: F945105P</p> <p>Original Lode Service part</p>  <p>7" Programmable Control Unit front panel</p>	<p>Drive roller for 938900 938901 939910 939906</p> <p>Partnumber: 938400</p> <p>Original Lode Service part</p>  <p>Drive roller for 938900 938901 939910 939906</p>	<p>Fan motor compartment</p> <p>Partnumber: 939708</p> <p>Original Lode Service part</p>  <p>Fan motor compartment</p>
--	---	--	---	--



<p>Frequency controller</p> <p>Partnumber: 250002</p> <p>Original Lode Service part</p>  <p>Frequency controller</p>	<p>Inclination motor</p> <p>Partnumber: 2200342</p> <p>Original Lode Service part</p>  <p>Inclination motor</p>	<p>Inclination sensor</p> <p>Partnumber: 9457433</p> <p>Original Lode Service part</p>  <p>Inclination sensor</p>	<p>Main switch with fuse</p> <p>Partnumber: 217201</p> <p>Lode Service Part</p>  <p>Main switch with fuse</p>	<p>Motor assembly</p> <p>Partnumber: 938440</p> <p>Original Lode Service part</p>  <p>Motor assembly</p>
<p>Motor Belt</p> <p>Partnumber: 932983</p> <p>Original Lode Service part</p>  <p>Motor Belt</p>	<p>Reverse roller 938900 938901 939910 939906</p> <p>Partnumber: 938410</p> <p>Original Lode Service part</p>  <p>Reverse roller 938900 938901 939910 939906</p>	<p>Treadmill Control Module 2 (TCM2)</p> <p>Partnumber: F945705</p> <p>Original Lode Service Part</p>  <p>Treadmill Control Module 2 (TCM2)</p>	<p>Walking belt 938900 938901 939910 939906</p> <p>Partnumber: 938241</p> <p>Original Lode Service part</p>  <p>Walking belt 938900 938901 939910 939906</p>	<p>Walking Deck 938900 938901 939910 939906</p> <p>Partnumber: 938221</p> <p>Original Lode Service part</p>  <p>Walking Deck 938900 938901 939910 939906</p>



15 - Specifications

Workload

Maximum speed	20 km/h	12.4 mph
Minimum operational speed	0.5 km/h	0.3 mph
Speed adjustment steps	0.1 km/h	0.1 mph
Positive elevation	25 %	
Elevation adjustment steps	0.5 %	

Comfort

Height handrail	90 cm	
Allowed user weight	225 kg	496 lbs
Length handrail	85 cm	

User Interface

Terminal operation mode	✓
-------------------------	---

Connectivity

Lode ERM interface protocol	✓
"HP Cosmos" interface protocol	✓
"Woodway" interface protocol	✓
"Trackmaster" interface protocol	✓
"GE" interface protocol	✓
USB connector	✓
RS232 in connector	✓

Dimensions

Walking surface length	150 cm	59.1 inch
Walking surface width	50 cm	19.7 inch
Step up height	17 cm	6.7 inch
Product length (cm)	215 cm	84.6 inch
Product width (cm)	80 cm	31.5 inch
Product height	130 cm	51.2 inch
Product weight	149 kg	328.5 lbs

Power requirements

Power in Horse Power	5.5 HP
Maximum rated power input	2500 VA
115 V AC 50/60 Hz (2 phases)	✓
230 V AC 50/60 Hz	✓
Maximum motor power	1.8 kW

Order info

Standards & Safety

IEC 60601-1:2005	✓
ISO 13485:2016 compliant	✓
ISO 9001:2015 compliant	✓
Standard emergency lanyard	✓

Certification

CE class Im according to MDD93/42/EEC	✓
CB according to IEC60950-1	✓

Environmental conditions

Minimum operational temperature	14 °C
Maximum operational temperature	40 °C
Minimum operational air pressure	70 kPa
Maximum operational air pressure	106 kPa
Minimum operational non-condensing humidity	30 %
Maximum operational non-condensing humidity	90 %
Minimum storage & transport temperature	-25 °C
Maximum storage & transport temperature	70 °C
Minimum air pressure storage & transport	50 kPa
Maximum storage & transport air pressure	106 kPa
Min. humidity storage & transport	10 %
Max. humidity storage & transport	95 %

Partnumber 939910